
**Amphetamines by Supported Liquid Extraction (AMPSLE)
for Analysis by LC-MS/MS**

1.0 Purpose - This procedure specifies the required elements for the extraction and quantitation of amphetamines using supported liquid extraction (SLE) for LC-MS/MS analysis.

2.0 Scope – This procedure applies to Toxicology in the Raleigh, Triad, and Western locations of the State Crime Laboratory.

3.0 Definitions – See Toxicology Definitions List

4.0 Equipment, Materials and Reagents

4.1 Equipment

- Centrifuge
- Mechanical Pipettes
- Class A Volumetric flasks
- TurboVap or equivalent evaporator

4.2 Materials

- Test tubes (16 x 125, 16 x 150 mm) with caps
- ISOLUTE SLE+ columns
- Glass Conical test tubes
- Vortexer
- Pipette tips
- LC vials with pre-slit septa caps

4.3 Reagents

- Ammonium Hydroxide
- Dichloromethane
- Isopropanol
- Hydrochloric Acid
- Methyl tert-butyl ether (MTBE)
- Methanol, HPLC grade or higher
- Water, HPLC grade or higher
- Negative Blood

4.4 Primary Reference Standards

• Amphetamine
• Methamphetamine
• 3,4 – Methylenedioxyamphetamine (MDMA)
• 3,4 – Methylenedioxyamphetamine (MDA)
• Mepivacaine

4.5 Prepared Reagents – Refer to [Toxicology Solution Prep Guidelines](#) for instructions on how to prepare the reagents required by this procedure.

4.6 Prepared Standards –Standards may be prepared by the Forensic Scientist in any amount provided that the component ratios are kept constant. If using 100 µg/mL reference standard solutions in place of a 1.0 mg/mL solution, multiply the volume pipetted by 10.

4.6.1 Mepivacaine Internal Standard Solution

4.6.1.1 Prepare the stock solution in methanol in a 50 ml volumetric flask. The final concentration of the mepivacaine will be 20 µg/mL.

4.6.1.1.1 Pipette 1 mL of the 1.0 mg/ml Mepivacaine solution into the volumetric flask.

4.6.1.1.2 QS to volume with methanol.

4.6.1.2 Lot number: Eight-digit format year/month/day

4.6.1.2.1 Example: 20230713

4.6.1.3 Expiration: One year.

4.6.1.4 Storage: Freezer.

4.6.1.5 QC check: Successful control check.

4.6.2 Amphetamines Internal Standard Working Solution (AISW)

4.6.2.1 Prepare the stock solution in methanol in a 5 ml volumetric flask. The final concentration of the mepivacaine will be 500 ng/mL.

4.6.2.1.1 Pipette 125 µL of Mepivacaine Internal Standard Solution (**4.6.1**) into a volumetric flask.

4.6.2.1.2 QS to volume with methanol.

4.6.2.2 Expiration: Prepare Daily.

4.6.2.3 QC check: N/A

4.6.3 Amphetamines Stock Calibration Solution (ASC)

4.6.3.1 Prepare the stock solution in methanol in a 50 ml volumetric flask. The final concentration of the amphetamines will be 10 µg/mL of Amphetamine, Methamphetamine, MDA, and MDMA.

4.6.3.1.1 Pipette 500 µL of each 1.0 mg/mL reference standard analyte solution into the volumetric flask.

4.6.3.1.2 QS to volume with methanol.

4.6.3.2 Lot number: Eight-digit format year/month/day + ASC

4.6.3.2.1 Example: 20230713ASC

4.6.3.3 Expiration: One year.

4.6.3.4 Storage: Freezer.

4.6.3.5 QC check: Successful calibration (see 5.3).

4.6.4 Amphetamines Calibration Working Solution (ACW)

4.6.4.1 Dilute 50 µL of the ASC solution with 950 µL of Methanol.

4.6.4.2 Expiration: Prepare Daily.

4.6.4.3 QC check: N/A

4.6.5 Amphetamines Stock Verification Solution (AV)

4.6.5.1 Prepare the stock solution in methanol in a 50 ml volumetric flask. The final concentration of the amphetamines will be 2.5 µg/mL.

4.6.5.2 The Verification Stock solution shall be prepared using standards from different manufacturers or different lot numbers from the ones used to prepare the calibration solution.

4.6.5.2.1 Pipette 125 µL of the 1.0 mg/mL reference standard analyte solution into the volumetric flask.

4.6.5.2.2 QS to volume with methanol.

4.6.5.3 Lot number: Eight-digit format year/month/day + AV

4.6.5.3.1 Example: 20230713AV

4.6.5.4 Expiration: One year.

4.6.5.5 Storage: Freezer.

4.6.5.6 QC check: Successful control check.

5.0 Procedure

5.1 Allow all solutions and samples to be analyzed to equilibrate to room temperature.

5.2 Prepare a labeled collection tube for each sample to be extracted and place in the appropriate location in the manifold.

5.2.1 Pipette 100 µL of 50 mM HCl in Methanol into each labeled collection tube.

5.3 Calibration and Control Sample Preparation

5.3.1 Calibrator Preparation

5.3.1.1 Pipette the following volumes of negative blood, the ASC, and the ACW into the appropriately labeled test tubes.

Amount of Neg. Blood (μL)	Amount of ASC (μL)	Amount of ACW (μL)	Final Concentration (ng/mL)
480		20	20
450		50	50
425		75	75
400		100	100
490	10		200
480	20		400
470	30		600

5.3.2 Positive Control Preparation

5.3.2.1 Pipette the following volumes of the AV and negative blood into an appropriately labeled test tube.

Amount of Neg. Blood (μL)	Amount of AV (μL)	Final Concentration (ng/mL)
480	20	60
340	160	480

5.2.3 Negative Control Preparation

5.2.3.1 Add 0.5 mL of negative blood to the labeled test tube.

5.2.4 An extraction batch will include at least two negative and two positive controls. Case specimens shall be bracketed by one of each.

5.2.5 Control samples shall make up at least 10 % of an extraction batch.

5.3 Maintenance

5.3.1 Add water to the TurboVap if needed.

5.4 Sampling

5.4.1 Pipette 0.5 mL of case specimens to be analyzed into a labeled and capped test tube.

5.4.1.1 Ensure that all body fluids are homogenous.

5.4.1.2 If a homogenous sample cannot be obtained, a notation shall be made in the worksheet detailing the condition of the sample and its handling.

5.5 Extraction Procedure

5.5.1 Add 100 μL of AISW to each calibrator, control, and case specimen to be analyzed.

5.5.2 Add 500 μL of 0.1 % ammonium hydroxide (aq) to each sample and vortex.

-
- 5.5.3 Centrifuge samples for five minutes.
 - 5.5.4 Pipette 750 µL of sample onto labeled ISOLUTE SLE+ column and allow to absorb onto the column.
 - 5.5.5 Add 2.5 mL 95:5 Dichloromethane/Isopropanol to each column and allow to elute.
 - 5.5.6 Add 2.5 mL MTBE to each column and allow to elute into the same tube as 5.5.5.
 - 5.5.7 Repeat 5.5.6.
 - 5.5.8 Evaporate to dryness using a TurboVap at 40 °C, 8-10 psi.
 - 5.5.9 Add 100 µL of methanol and vortex.
 - 5.5.10 Add 400 µL of Mobile Phase A and vortex.
 - 5.5.11 Centrifuge for 5 minutes.
 - 5.5.12 Transfer reconstituted specimens to labeled LC vials and cap.
 - 5.5.13 Analyze samples on a LC-MS/MS as specified in the [Toxicology Liquid Chromatography-Tandem Mass Spectrometry \(LC-MS/MS\) procedure](#).

5.6 Data Processing and Calibration/Control Acceptance Criteria

- 5.6.1 Process the run using the AMPSLE method.
 - 5.6.1.1 Ensure that the boxes are checked to update the retention times and ion ratios of the analytes and their internal standards.
- 5.6.2 Examine the calibration samples for outliers. See 5.6.5.1 and 5.6.5.2.
- 5.6.3 If manual integration of a compound is needed, the chromatogram showing the integration prior to manual integration shall be printed and included with the data as well. The reason for the manual integration shall be documented on the chromatogram.
- 5.6.4 Save the data with the name of the procedure and the extraction date added to the end.
 - 5.6.4.1 Example: AMPSLE20250428
- 5.6.5 **Calibration Curve Acceptance Criteria**
 - 5.6.5.1 Evaluate the curve by back calculating the calibrator concentrations against the curve. Values of +/- 25 % from the target concentration are acceptable for the lowest calibrator. All other calibrators shall be within 20 % of the target concentration.
 - 5.6.5.2 The qualifier ion ratios for an analyte/internal standard in the calibrators shall be within +/- 20 % of the analyte ion ratios determined by the average of the calibration sample ion ratios.

-
- 5.6.5.3** A maximum of two calibration points may be dropped from the curve with cause (e.g., statistical outlier, laboratory accident, fails **5.6.5.1**, **5.6.5.2**, etc.).
- 5.6.5.3.1** If the low calibration point is dropped, this will change the LOQ and may require a repeat analysis for case specimens whose quantitation is between the lowest extracted calibrator and the new LOQ.
- 5.6.5.3.2** If the high calibration point is dropped, this will change the ULOQ and may require a repeat analysis for case specimens whose quantitation is between the highest extracted calibrator and the new ULOQ.
- 5.6.5.3.3** If the two lowest or two highest calibration points for an analyte are dropped, the run fails for that analyte and the extraction must be repeated.
- 5.6.5.4** The calibration curves for each analyte shall show a coefficient of determination (r^2) of 0.985 or greater.
- 5.6.5.5** The absolute value of the x-intercept shall be less than the lowest extracted calibrator.
- 5.6.5.6** If a calibration curve for an analyte fails to meet the criteria in **5.6.5.4** or **5.6.5.5**, the extraction shall be repeated.

5.6.6 Quality Control Acceptance Criteria

- 5.6.6.1** Each analyte in a positive control shall give a quantitation within +/- 20 % of the expected concentration.
- 5.6.6.2** The qualifier ion ratios for an analyte/internal standard in the controls shall be within +/- 20 % of the analyte ion ratios determined by the average of the calibration sample ion ratios.
- 5.6.6.3** The negative control fails for an analyte if there is an integrated peak for both transitions at the expected retention time that is greater than 20 % of the lowest calibrator area of that analyte and meets the requirement in **5.6.6.2**.
- 5.6.7** The failure of an analyte to meet the criteria in **5.6.5** and **5.6.6** does not invalidate the acceptability of another analyte.
- 5.6.8** Create a data packet for the run, including the following quality control data:
- Summary page with FA workstation reference
 - Completed extraction worksheet.
 - LC-MS/MS sequence list
 - Approved LC-MS/MS system check
 - Experiment, Method, and Calibration Report
 - Quantitation reports of all calibrators and controls

5.6.9 The Quality Control data packet will be named beginning with “AMPSLE” (capitalization optional), followed by eight-digit format year/month/day ending with the instrument name. A suffix may be added to differentiate multiple runs.

5.6.9.1 Example: AMPSLE20251004-XXX-R-X

5.6.10 All quality control data packets shall be administratively and technically reviewed prior to use of the associated case data for reporting.

5.6.11 The reviewed data packet shall be uploaded into the workstation managed files. The review and approval will be indicated by signing the summary page prior to uploading to FA.

5.6.12 Control Charting

5.6.12.1 Complete the Toxicology Control Chart Form.

5.7 Sample Acceptance Criteria

5.7.1 Analyte and Internal Standard Identification Criteria

5.7.1.1 The qualifier ion ratios shall be within +/- 20 % of the target value.

5.7.1.2 The retention time shall not differ from the target value by more than 2.0 %.

5.7.2 Quantitative Acceptance Criteria

5.7.2.1 The internal standard area must be within 50 % - 200 % of the average internal standard area of the controls.

5.7.2.1.1 If the internal standard area is less than 50 %, the quantitative results may be used if the signal to noise (S/N) calculated by the instrument software is greater than 10:1. If not, the case shall be re-analyzed to confirm the quantitation, sample volume permitting.

5.7.2.1.2 If the internal standard area is greater than 200 %, the case shall be re-analyzed to confirm the quantitation, sample volume permitting.

5.7.2.1.3 If the internal standard area fails to meet acceptance criteria and there is insufficient volume remaining, the data may be evaluated qualitatively with documented approval from the Toxicology Technical Leader.

5.7.2.2 The signal to noise (S/N) of each analyte to be quantitated shall be greater than 10:1.

5.7.2.3 The quantitation result shall be equal to or greater than the LOQ of each analyte to be reported.

5.7.3 Qualitative Acceptance Criteria

-
- 5.7.3.1** If there is insufficient sample volume remaining, the data may be reported qualitatively only if the following acceptance criteria are met:
- 5.7.3.1.1** Calibration curve is acceptable.
 - 5.7.3.1.2** Quality control acceptance criteria in **5.6.6.2** and **5.6.6.3**.
 - 5.7.3.1.3** Case sample acceptance criteria in **5.7.1** and **5.7.2**.
- 5.7.4** If an analyte meets the criteria in **5.7.1** and the corresponding internal standard fails to meet the same criteria; the case sample will be re-analyzed, sample volume permitting. If there is insufficient volume remaining, the data may be evaluated qualitatively with documented approval from the Toxicology Technical Leader.
- 5.7.5** If an analyte meets the criteria listed in **5.7.2**, but fails to meet the criteria in **5.7.1**, the case sample will be re-analyzed. The data may not be used qualitatively for that analyte.
- 5.7.6** If an analyte meets the criteria in **5.7.1** and **5.7.2** but failed to meet the acceptance criteria listed in **5.6.6.1**, the case shall be re-analyzed, sample volume permitting.
- 5.7.6.1** If there is insufficient volume remaining the data may be reported qualitatively, provided all acceptance criteria in **5.7.3** are met.
- 5.8 Calculations**
- 5.8.1** Reported Measurement uncertainty calculation: see **5.10.6**.
 - 5.8.2** Percent Difference Calculation:
$$\frac{|(\text{standard retention time} - \text{analyte retention time})|}{(\text{standard retention time})} * 100$$
- 5.9 Uncertainty of Measurement**
- 5.9.1** The current process uncertainty for each analyte quantitated by this procedure is located in the [Toxicology Reporting Index](#).
 - 5.9.2** The uncertainty of measurement shall be reported in the same units as the compound quantitation.
- 5.10 Reporting**
- 5.10.1** Amphetamines identified by LC-MS/MS analysis shall have a positive indication from an initial analysis of another aliquot. Refer to the [Drug Toxicology Reporting](#) procedure for reporting of amphetamines.
 - 5.10.2** All quantitative results will be truncated and reported to two significant figures.
 - 5.10.3** All calculated measurement uncertainties will be rounded and reported to the same level of significance as the quantitative results.
 - 5.10.4** If a case sample is re-extracted, the average of all acceptable quantitations will be reported.

-
- 5.10.5** The full analyte result is used to calculate the associated measurement uncertainty (with coverage factor $k=2$, 95.45% coverage probability).
- 5.10.6** Example Reporting Statement and Uncertainty determination.
- 5.10.6.1** Methamphetamine full quantitation = 50.25 ng/mL
- 5.10.6.2** Multiply the full quantitation by the current process uncertainty (see 5.9.1).
- 5.10.6.2.1** $50.25 \text{ ng/mL} * 0.50 (50\%) = 25.125 \text{ ng/mL}$.
- 5.10.6.3** The Methamphetamine result is truncated to 50 ng/mL, and the associated uncertainty is rounded to 25 ng/mL for reporting.
- 5.10.6.4** Report: Methamphetamine – 50 ng/ml +/- 25 ng/mL at a coverage probability of 95.45%
- 5.10.7** Case Samples with concentrations that exceed the upper level of the calibration curve will be reanalyzed with a lower sample volume or after dilution with the proper matrix to bring it within the calibration range.

5.11 Record the following in the case record:

- Quantitation report of the sample

6.0 Limitations

6.1 Stability – Processed samples were determined to be stable for 52 hours after recon (**5.5.9**)

6.2 Dilution Integrity – 1:2 and 1:5 dilutions were validated for this method

6.3 Interferences – No interferences were determined during validation

7.0 Safety

7.1 Refer to the [Laboratory Safety Manual](#).

8.0 References

Biotage Application Note AN875 – Extraction of a Drugs of Abuse Pane from Whole Blood Using Isolute SLE+ Prior to UPLC-MS/MS Analysis

AMPSLE Method Validation

9.0 Records

- Quality Control data packet
- Case record
- Toxicology Control Chart Form

10.0 Attachments

-
- **Appendix 1** – Tune Settings
 - **Appendix 2** – Chemical Hygiene and Safety Precautions for Extremely Hazardous and Particularly Hazardous Substances

Revision History		
Effective Date	Version Number	Reason
08/01/2025	1	Original Document

Appendix 1 Tune Settings

1.0 Analyser settings - Shall be updated with each new Calibration and Resolution.

1.1 See [Toxicology Liquid Chromatography Tandem Mass Spectrometry \(LC-MS/MS\)](#) for Calibration – Mass Scale and Resolution requirements

1.2 The AMPSLE procedure utilizes a 0.75 Da Calibration and Resolution

1.3 The Calibration and Resolution file shall be linked to the tune file

The screenshot displays the Waters Xevo TQ-S cronos MS Detector software interface. The window title is "Waters Xevo TQ-S cronos MS Detector - D:\20250512.PRO\ACQUDB\AMPSLE20250117.ipr". The interface includes a menu bar (File, View, Ion Mode, Calibration, Gas, Vacuum, Ramps, Setup, Acquire, Help) and a toolbar. The main area is divided into two sections: a left-hand control panel and a right-hand data display area.

Control Panel (Left):

- Source Voltages:** Capillary (kV) set to 0.09, 1.00; Cone (V) set to -1, 50.
- Source Temperatures:** Desolvation Temp (°C) set to 59, 400.
- Source Gas Flow:** Desolvation (L/Hr) set to 16, 800; Cone (L/hr) set to 6, 20.
- Analyser:** LM Resolution 1 (10.2), HM Resolution 1 (15.0), Ion Energy 1 (-0.5); LM Resolution 2 (8.1), HM Resolution 2 (13.4), Ion Energy 2 (0.5); Collision Energy MS (V) (3).

Data Display Area (Right):

Function	Set	Mass	Span	Gain
<input checked="" type="checkbox"/> 1 MS Scan	56	304	10	10
<input checked="" type="checkbox"/> 2 MS Scan	219	318	10	10
<input checked="" type="checkbox"/> 3 MS Scan	502	290	10	10
<input type="checkbox"/> 4 MS1 Scan	614	922.35	5	9

Below the table are three mass spectra plots for masses 304.0, 318.0, and 290.0. Each plot shows a single sharp peak at the respective mass value with a scale of 1.00e0 x10. The x-axis for the 304.0 plot ranges from 300.0 to 305.0, for the 318.0 plot from 315.0 to 320.0, and for the 290.0 plot from 5.0 to 290.0.


At the bottom of the interface, there are status indicators: "Ready", "Vacuum Ok", and "Standby".

The screenshot displays the Waters Xevo TQ-S cronos MS Detector software interface. The window title is "Waters Xevo TQ-S cronos MS Detector - D:\20250512.PRO\ACQUDB\AMPSLE20250117.ipr". The menu bar includes File, View, Ion Mode, Calibration, Gas, Vacuum, Ramps, Setup, Acquire, and Help. The interface is divided into several sections:

- Source Parameters:** Source Temp (°C) is set to 147, with a target of 120. Ion Guide 2 Offset (V) is 0.5. Collision Cell Lenses: Entrance is 1.7, Exit is -0.7. DETECTOR Gain is 1.00. A "Reset to Defaults" button is present.
- Scan Table:** A table with columns for Function, Set, Mass, Span, and Gain.

Function	Set	Mass	Span	Gain
<input checked="" type="checkbox"/> 1 MS Scan	56	304	10	10
<input checked="" type="checkbox"/> 2 MS Scan	219	318	10	10
<input checked="" type="checkbox"/> 3 MS Scan	502	290	10	10
<input type="checkbox"/> 4 MS1 Scan	614	922.35	5	9
- Mass Spectra:** Three panels showing mass spectra for m/z 304.0, 318.0, and 290.0. Each panel has a y-axis scale of 1.00e0 x10 and an x-axis from 300.0 to 305.0, 315.0 to 320.0, and 5.0 to 29.0 respectively.
- Status Bar:** Shows "Ready", "Vacuum Ok", and "Standby" indicators.

Appendix 2:
 Chemical Hygiene and Safety Precautions for Extremely Hazardous and Particularly Hazardous Substances

Methylene Chloride/Dichloromethane							
DANGER: PARTICULARLY HAZARDOUS SUBSTANCE							
	<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <td style="background-color: #0056b3; color: white;">HEALTH</td> <td style="text-align: right; color: white;">2</td> </tr> <tr> <td style="background-color: #ff0000; color: white;">FLAMMABILITY</td> <td style="text-align: right; color: white;">1</td> </tr> <tr> <td style="background-color: #ffff00;">REACTIVITY</td> <td style="text-align: right;">1</td> </tr> </table>	HEALTH	2	FLAMMABILITY	1	REACTIVITY	1
HEALTH	2						
FLAMMABILITY	1						
REACTIVITY	1						
Detection of Release	Clear colorless liquid. Ether like odor						
Signs/Symptoms of Exposure	Serious eye irritation; skin irritation; may cause drowsiness or dizziness.						
PEL	ACGIH (TLV) – 50 ppm; OSHA Specifically Regulated Chemicals/Carcinogens – (PEL) 25 ppm						
Associated Hazards	Serious eye and skin irritation; suspected of causing cancer						
Controls	Use under fume hood. Avoid contact with skin, eyes and clothing. Wash hands before breaks and immediately after handling the product. Use eye protection. Handle with gloves. Wear lab coat. Gloves: Fluorinated rubber (break through time = 148 minutes)						
Safe handling, storage, disposal	Avoid contact with skin and eyes. Avoid inhalation of vapor or mist. Keep in a tightly closed container. Containers which are opened must be carefully resealed and kept upright to prevent leakage. Dispose of in Hazardous Chemical Waste.						
Emergency Procedures (2.2)(4.1)(6)	<p>Eye Contact: Rinse thoroughly with plenty of water for at least 15 minutes and consult a physician.</p> <p>Inhalation Exposure: If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician.</p> <p>Ingestion: Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.</p> <p>Skin Contact: Wash off with soap and plenty of water. Consult a physician.</p> <p>Spills: Avoid breathing vapors, mist or gas. Ensure adequate ventilation. Evacuate personnel to safe areas. Small contained spill: wearing appropriate PPE, collect with absorbent material, and place in container. Dispose in Hazardous Chemical Waste. Large spills: Evacuate area and call 911 (Haz Mat).</p>						