

## **Technical Procedure for the Examination of Primer Gunshot Residue Evidence**

- 1.0 Purpose** – This technical procedure shall be followed for the examination of primer gunshot residue evidence.
- 2.0 Scope** – This procedure applies to gunshot residue (GSR) kits, and adhesive lifts collected from inanimate objects.
- 3.0 Definitions** – N/A
- 4.0 Equipment, Materials, and Reagents**
- 4.1 Equipment**
- Scanning Electron Microscope/Energy Dispersive X-Ray Detector (SEM/EDS)
- 4.2 Materials**
- Adhesive lifts
- 4.3 Reagents**
- 10% bleach solution
- 5.0 Procedure**
- 5.1** Antimony (Sb), barium (Ba), and lead (Pb) are the three most common elements associated with primer gunshot residue.
- 5.2** Each GSR collection kit shall contain adhesive lifts used to collect evidence from the hands of an individual or inanimate objects (e.g., clothing, vehicles, or other surfaces). Any additional collection materials present (e.g., swabs) will not be examined.
- 5.3 Analytical Approach for GSR kit examination**
- 5.3.1** Review the request for analysis.
- 5.3.2** Open kit. Label each piece of evidence within the kit with the item number, date, and the Forensic Scientist's initials.
- 5.3.3** If kit paperwork (e.g., GSR Analysis Information Form) is present in the kit, label the top of the form with the Laboratory case number, Laboratory item number, Forensic Scientist's initials, and the date.
- 5.3.4** Scan the labeled paperwork into Forensic Advantage (FA). If this paperwork is not present in the kit, note this in the GSR case record worksheet.
- 5.3.5** Determine if the GSR kit is suitable for examination. GSR kits from the hands of an individual that meet one or more of the following criteria shall not be examined.



## **5.6 Guidelines for primer gunshot residue results statements**

- 5.6.1** All results shall be based on the Forensic Scientist's knowledge and experience and the case being examined. The Forensic Scientist has the responsibility to conduct analyses deemed necessary and appropriate and report the results based on his/her findings and experience.
- 5.6.2** The following is a list of guidelines for reporting results of primer gunshot residue analysis:
- 5.6.2.1** Examination of the adhesive lift(s) collected from item(s) # \_\_\_\_\_ revealed the presence of particle(s) characteristic to primer gunshot residue. Characteristic particle(s) on an item could have originated from being in close proximity to a firearm when it was discharged or coming into contact with another item with primer gunshot residue on it.
  - 5.6.2.2** Examination of the adhesive lifts collected from Item(s) # \_\_\_\_\_ did not identify particles characteristic of primer gunshot residue.
  - 5.6.2.3** Examination of the gunshot residue kit revealed the presence of particle(s) characteristic to primer gunshot residue. Characteristic particle(s) could have originated from the discharge of a firearm, the handling of a discharged firearm, being in close proximity to a firearm when it was discharged or from contacting an object with characteristic primer gunshot residue particles on it.
  - 5.6.2.4** Examination of the gunshot residue kit did not identify particles characteristic of primer gunshot residue. There is no indication the subject came into contact with particles characteristic of primer gunshot residue or was in the vicinity of a firearm during discharge. The absence of characteristic particles of primer gunshot residue on hands does not conclude that a subject did not fire a gun. It is possible the subject fired a gun, but the characteristic particles of primer gunshot residue were not deposited, were removed by activity, or were not detected.
  - 5.6.2.5** Examination of the case details revealed that the subject had washed his/her hands prior to the gunshot residue kit being administered. Washing one's hands removes primer gunshot residue. Therefore, the evidence was not examined.
  - 5.6.2.6** Examination of the case details revealed that a time greater than (4) hours had elapsed from the time the weapon was discharged to the time the evidence was collected. The concentration of primer gunshot residue significantly declines after approximately (4) hours. Therefore, the evidence was not examined.
  - 5.6.2.7** This item was not examined due to an improper collection technique. Proper collection for characteristic primer gunshot residue particles requires the use of adhesive lifts.

**5.6.2.8** Examination of item(s) # \_\_\_\_\_ revealed the presence of characteristic particles of primer gunshot residue; therefore, this item was not analyzed.

**5.6.2.8.1** This statement shall be used when characteristic particle(s) of primer gunshot residue are identified on an item collected from the same subject or inanimate object.

**5.6.2.9** Characteristic particles of primer gunshot residue can be present on a subject that sustained a gunshot wound. Therefore, the GSR kit/clothing was not examined.

**5.6.2.10** Spheroidal \_\_\_\_\_ (indicate elemental composition) particles were found on the adhesive lifts collected from \_\_\_\_\_ (Item #); therefore the \_\_\_\_\_ (Item #) was analyzed for comparison. The elemental composition of the particles on the adhesive lifts collected from \_\_\_\_\_ (Item #) were consistent/not consistent with the elemental composition of the particles from \_\_\_\_\_ (Item #). These particles, however, could be found in [other] ammunitions or [other] sources.

**5.7** **Calibrations** – No additional calibrations or performance checks are required. See associated technical procedures for instrumental performance checks.

**5.8** **Maintenance** – No additional maintenance is required. See associated technical procedures for instrumental maintenance.

**5.9** **Sampling**

**5.9.1** No sampling is performed. When sample selection occurs, it shall be based on the Forensic Scientist’s training and experience.

**5.10** **Gunshot Residue Examination Criteria**

**5.10.1** Items submitted for primer gunshot residue examination shall be analyzed based on the following criteria:

**5.10.1.1** When an agency requests gunshot residue examination for both a gunshot residue kit and clothing from the same subject, the gunshot residue kit shall be examined first.

**5.10.1.1.1** If examination of the gunshot residue collection kit reveals particle(s) characteristic of primer gunshot residue, the clothing shall not be examined.

**5.10.1.1.2** If examination of the gunshot residue collection kit does not reveal the presence of particle(s) characteristic of primer gunshot residue, the clothing shall be examined.

**5.10.1.2** Gunshot residue kits and clothing from a subject that sustained a gunshot wound shall not be examined in the Trace Evidence Section for characteristic particles of primer gunshot residue.

**5.10.1.3** Inanimate objects on which primer gunshot residue is expected will not be analyzed. Examples of items that will not be analyzed for gunshot residue include, but are not limited to, firearms and inanimate objects with known bullet holes.

**5.11 Calculations** – N/A

**5.12 Uncertainty of Measurement** – N/A

## **6.0 Limitations**

**6.1** The presence of characteristic particles of primer gunshot residue on a subject's hands or clothing does not conclude that the subject fired a gun. Characteristic particles could have originated from the discharge of a firearm, the handling of a discharged firearm, being in close proximity to a firearm when it was discharged, or from contacting an object with characteristic primer gunshot residue particles on it.

**6.2** The absence of characteristic particles of primer gunshot residue on a subject's hands or clothing does not conclude that a subject did not fire a gun. It is possible the subject fired a gun, but the characteristic particles of primer gunshot residue were not deposited, were removed by activity, or were not detected.

**6.3** Inanimate objects packaged together shall be treated as one item. Primer GSR particles can be transferred between objects if packaged together.

## **7.0 Safety**

**7.1** Items of clothing or other inanimate objects may have blood or other body fluids present. Use proper protective equipment when dealing with items that may contain biohazard material.

**7.2** X-Rays from the EDS detector of the scanning electron microscope are considered to be a minimal health risk. Dosimeter badges are located near the systems to monitor radiation levels.

## **8.0 References**

ASTM Standard E1588, "Standard Practice for Gunshot Residue Analysis by Scanning Electron Microscopy/ Energy—Dispersive X-ray Spectrometry." ASTM International, West Conshohocken, PA, [www.astm.org](http://www.astm.org).

ASTM Standard E3309, "Standard Guide for Reporting of Forensic Primer Gunshot Residue (pGSR) Analysis by Scanning Electron Microscopy/Energy Dispersive X-Ray Spectrometry (SEM/EDS)." ASTM International, West Conshohocken, PA, [www.astm.org](http://www.astm.org).

Andrasko, J. "Detection of Gunshot Residue on Hands by Scanning Electron Microscopy." *Journal of Forensic Sciences* 22.2 (1977): 279-287.

Basu, Samarendra. "Formations of Gunshot Residues." *Journal of Forensic Sciences* 27.1 (1982): 72-91.

DeGaetano, Douglas, et al. "A Comparison of Three Techniques Developed for Sampling and Analysis of Gunshot Residue by Scanning Electron Microscopy/Energy Dispersive X-Ray Analysis (SEM/EDX)." *Journal of Forensic Sciences* 37.1 (1992): 281-300.

Garofano, L., et al. "Gunshot Residue Further Studies on particles of Environmental and Occupational Origin." *Forensic Science International* 103 (1999): 1-21.

Kilty, J.W. "Activity After Shooting and Its Effect on the Retention of Primer Residue." *Journal of Forensic Sciences* 20.2 (1975): 219-230.

Molina, Kimberly D., et al. "Gunshot Residue Testing in Suicides Part I: Analysis by Scanning Electron Microscopy with Energy-Dispersive X-Ray." *The American Journal of Forensic Medicine & Pathology* 28.3 (2007): 187-190.

Nesbitt, R.S., et al. "Detection of Gunshot Residue by Use of the Scanning Electron Microscope." *Journal of Forensic Sciences* 21.3 (1976): 595-610.

Romolo, Francesco Saverio, et al. "Identification of Gunshot Residue: A Critical Review." *Forensic Science International* 119 (2001): 195-211.

Schwoeble, A.J. & David Exline. *Current Methods in Forensic Gunshot Residue Analysis*. Boca Raton: CRC Press, 2000.

SWGGSR. "Guide for Primer Gunshot Residue Analysis by Scanning Electron Microscopy/Energy Dispersive X-Ray Spectrometry 11-29-11." [www.swggsr.org](http://www.swggsr.org)

Wallace, Jim. "Chemical Aspects of Firearms Ammunition." *AFTE Journal* 22.4 (1990): 364-389.

Wolten, G.M., et al. "Final Report on Particle Analysis for Gunshot Residue Detection." The Aerospace Corporation, ATR-77 (7915)-3, 1977.

Wright, Diana M., & Trimpe, Michael A. "Summary of the FBI Laboratory's Gunshot Residue Symposium, May 31-June 5, 2005." *Forensic Science Communications* 8.3 (2006).

Zeichner, Arie. "More on the Uniqueness of Gunshot Residue (GSR) Particles." *Journal of Forensic Sciences* 42.6 (1997): 1027-1028.

**9.0 Records – N/A**

**10.0 Attachments – N/A**

<b>Revision History</b>		
Effective Date	Version Number	Reason
07/07/2023	1	Original document created from the Technical Procedure for the Examination of Gunshot Residue Evidence