

Toxicology Unit

Effective Date: 10/05/2015

Pipette Number	Calibration Performance Check <i>Select either Fixed or Adjustable volume pipette</i>	QC Criteria Results	Comments/Maintenance
Date: _____	Fixed - Volume tested: _____ 1. _____ grams 2. _____ grams _____ Adjustable - Lowest volume setting tested: _____ 1. _____ grams 2. _____ grams Adjustable - Highest volume setting tested: _____ 1. _____ grams 2. _____ grams	All volumes within 3% for >0.01g? Yes _____ No _____ NA _____ All volumes within 10% for <0.01g Yes _____ No _____ NA _____	
Initials: _____			
Date: _____	Fixed - Volume tested: _____ 1. _____ grams 2. _____ grams _____ Adjustable - Lowest volume setting tested: _____ 1. _____ grams 2. _____ grams Adjustable - Highest volume setting tested: _____ 1. _____ grams 2. _____ grams	All volumes within 3% for >0.01g? Yes _____ No _____ NA _____ All volumes within 10% for <0.01g Yes _____ No _____ NA _____	
Initials: _____			