

QCS	QTOF STANDARDS	LOD (ng/ml)	Cutoff (ng/ml)	Notes
3	Aripiprazole	5	25	
3	Atomoxetine	40	50	
3	Chlordiazepoxide	5	10	
3	Chlorpheniramine	20	20	
3	Clonidine	10	10	
3	Clozapine	50	50	
3	Cocaethylene	2.5	10	
3	Flunitrazepam	5	5	
3	Isobutryl Fentanyl	2.5	2.5	
3	Loratadine	10	10	
3	Meclizine	40	50	
3	Meperidine	25	25	
3	Methylone	40	50	
3	Mirtazapine	20	25	
3	Noroxycodone	50	50	
3	Oxazepam	40	40	
3	Pentylone	40	50	
3	Promethazine	25	25	
3	Thioridazine	25	50	
3	Triazolam	5	5	
3	Venlafaxine	20	20	
3	Zolpidem	2.5	10	
4	3 4-Methylenedioxypyrovalerone (MDPV)	20	25	
4	7-Aminoclonazepam	5	10	
4	Alpha-Pyrrolidinovalerophenone (Alpha-PVP)	40	40	
4	Alprazolam	20	20	
4	Amoxapine	10	100	
4	Butylone	20	20	
4	Carisoprodol	25	400	
4	Desalkylflurazepam	30	30	
4	Desipramine	100	100	
4	Diclazepam	5	10	
4	Doxepin	20	50	
4	Etizolam	10	10	
4	Methoxetamine	5	5	
4	Metoclopramide	20	50	
4	Oxymorphone	40	40	
4	Procainamide	25	500	
4	Quetiapine	2.5	10	
4	Trazodone	25	25	
4	U-47700	40	40	
4	Verapamil	10	25	

[illegible]

[illegible]

[illegible]

[illegible]

--	--	--	--	--